



		RAMP LENGTH, CM	TRANSITION HEIGHT, CM			
<b>Personal</b>		70	5	9	14	18
	<b>Multipurpose</b>	87.5	6	12	18	23
		100	7	14	21	26
<b>Senior</b>	<b>EBL</b>	90	6	13	19	23
		120	8	17	25	31
		135	9	19	28	35
		150	10	21	31	39
		165	12	23	34	43
		200	14	28	42	52
		250	17	35	52	65
		300	21	42	62	78
		350	24	49	73	91
		185	13	26	38	48

Invacare recommends risk assessment is completed for use over 4 degrees

AS/NZ standards only apply to 4 degrees

\* The Ramp Calculator is indicative only to assist your judgement